

KONTED



C10CX

***LINEAR Color Doppler***

[www.konted.cn](http://www.konted.cn)



Frequency

Depth

Mode

Applications

7.5-10MHz

20-100mm

B Mode  
 B+M  
 PW Doppler  
 Color Doppler  
 Power Doppler

MSK  
 Breast  
 Thyroid  
 Superficial  
 Vascular  
 Pediatrics  
 Carotid  
 Small Parts

# Download our ultrasound software from App Store or Google Play



- ✓ Unlimited users
- ✓ Free software
- ✓ No hidden fees
- ✓ Support printing reports
- ✓ Chinese, English, Russian, Italian  
Spanish, Portuguese(Brazil)

## Software Controls

GN(gain)

D(depth)

ENH(enhancement)

DR(dynamic range)

F(frequency)

FocusPos

PRF

WF

Mode

8TGC

Biopsy

Annote

Page left/right up/down

## Measurement

B: Length, Area/Circumference,  
Angle, Trace, Distance

B+M: Heart Rate, Time, Distance

B+PW: Velocity, Heart Rate(2), S/D, Depth



Package includes:

- 1\* wireless probe
- 1\* wireless charger
- 1\* USB charger
- 1\* wrist strap
- 1\* plug
- 1\* portable aluminium alloy box

## Connection Type

*WiFi (802.11n/20MHz/5G)  
Probe built-in WiFi signal*

## Supporting Devices

*Tablet, Smart phone, PC, etc*

## Working system

*Android, IOS, Windows*

## Image save

*jpg, avi, Mp4 and DICOM format*

## Buttons

*3 programmable buttons*

## Lithium Battery

*4200 mAh. 3.8V Support wireless charging  
Working time: ≥3h, Standby time: ≥12h  
Battery life 5-8years*

## Dimensions

*158(h)x70(w)x30(d) mm*

Net weight 223g

Elements 192


Footprint L45

## Boot-up

*Platform dependent, generally <5 seconds*



# Offer great images to everyone





- \* **Training how to install the machine, training how to use the machine.**
- \* **Spare parts replacement and other software/hardware trouble shooting**
- \* **Offer OEM&ODM service for distributors**
- \* **Offer SDK depends on different models**
- \* **Free software and free upgrade the software**
- \* **18 Months Warranty, extendable up to 3 years**

**BEIJING KONTED MEDICAL TECHNOLOGY CO.,LTD**

**Tel:8610-60219113**

**Fax:8610-60219213**

**Service and Cooperation :info@konted.cn**

**Website:www.konted.cn**

**Facebook: Beijing Konted**

**Youtube: Medical Konted**

**Address: Room 111,1F, Building 3, No. 27 Yongwang Road,Daxing Biological Pharmaceutical Industry Base,Daxing District,Beijing 102629,China.**